

MONOFOKAL

HARMONY



Ciljna grupa

Monofokalno sočivo sa blago pojačanim dodatkom u donjem delu sočiva. Preporučuje se kategoriji ljudi koji osećaju napor, tenziju u očima kao i glavobolju usled dužeg rada sa digitalnim uređajima.

Dizajn 0.50 – starosna dob 20+
Dizajn 0.75 – starosna dob 30+

Tehnologija

OPTOTECH Mild Add

Način izrade

Unutrašnja izrada monofokalne površine.

Poručivanje dioptrija

1/4

Minimalna visina ugradnje

18 mm

Optimizacija

Specijalno tanjenje za sočivo u plus opsegu na osnovu oblika okvira i parametra za fitovanje sočiva u okvir. Ukoliko se sočivo poručuje bez opcije Optimizacija, topografija sočiva će biti decentrirana nazalno u odnosu na geometrijski centar za 2,5 mm.

Ocene korisnika

Vid – daljina



Vid – blizina



Spontano privikavanje



Vid – srednja daljina



Komfor pri korištenju




Korišćenje digitalnih uređaja



Oznake na sočivima

- gravure
- markeri



MATERIJALI

| | |
|-----------------|---|
| 1.50 CR-39 | 1 |
| 1.53 Trivex | 2 |
| 1.56 Mid. index | 3 |
| 1.59 Poly | 4 |
| 1.60 High Index | 5 |
| 1.60 Tribid | 6 |
| 1.67 High Index | 7 |
| 1.74 High Index | 8 |

DORADE

| | | | | | |
|-----------------------------|---|-----------------------------|---|--------------------------------------|---|
| Bez dorade | W | Transitions GEN 8 Sapphire | F | Nupolar Gray | S |
| Remuve 400 | H | Transitions GEN 8 Amethyst | M | Nupolar Gray Light | L |
| Remuve 420 | R | Transitions GEN 8 Amber | A | Nupolar Green | E |
| Remuve 420 Photomatic Brown | O | Transitions GEN 8 Emerald | E | Transitions XTRActive Polarised Gray | X |
| Remuve 420 Photomatic Gray | U | Transitions XTRActive Brown | C | Photomatic Brown | B |
| Transitions GEN 8 Brown | Z | Transitions XTRActive Gray | L | Photomatic Gray | I |
| Transitions GEN 8 Gray | N | Transitions XTRActive Green | K | DriveWear | D |
| Transitions GEN 8 Green | V | Nupolar Brown | B | | |