

PROGRESIV

POLARISMD



Ciljna grupa

Iskusni poznavaooci progresivnih sočiva. Starosna dob 45+.

Tehnologija

OPTOTECH MD
The Universal Vision

Način izrade

Unutrašnja izrada progresivne površine.

Poručivanje dioptrija

1/4

Progresivna dužina

15 mm, 17 mm

Minimalna visina ugradnje

17 mm, 19 mm

Optimizacija

Specijalno tanjenje za sočivo u plus opsegu na osnovu oblika okvira i parametara za fitovanje sočiva u okvir. Ukoliko se sočivo poručuje bez opcije Optimizacija, sočivo će biti decentrirano za 2,5 mm nazalno.

Ocene korisnika

Vid – daljina



Vid – blizina



Spontano privikavanje



Vid – srednja daljina



Komfor pri korištenju




Korištenje digitalnih uređaja



Oznake na sočivima

- gravure
- markeri



MATERIJALI

1.50 CR-39	1
1.53 Trivex	2
1.56 Mid. index	3
1.59 Poly	4
1.60 High Index	5
1.60 Tribrid	6
1.67 High Index	7
1.74 High Index	8

DORADE

Bez dorade	W	Transitions GEN 8 Sapphire	F	Nupolar Gray	S
Remove 400	H	Transitions GEN 8 Amethyst	M	Nupolar Gray Light	L
Remove 420	R	Transitions GEN 8 Amber	A	Nupolar Green	E
Remove 420 Photomatic Brown	O	Transitions GEN 8 Emerald	E	Transitions XTRActive Polarised Gray	X
Remove 420 Photomatic Gray	U	Transitions XTRActive Brown	C	Photomatic Brown	B
Transitions GEN 8 Brown	Z	Transitions XTRActive Gray	L	Photomatic Gray	I
Transitions GEN 8 Gray	N	Transitions XTRActive Green	K	DriveWear	D
Transitions GEN 8 Green	V	Nupolar Brown	B		